

Hiden Analytical Ltd. 420 Europa Boulevard Warrington WA5 7UN England

- T +44 [0] 1925 445 225
- F +44 [0] 1925 416 518
- E info@hiden.co.uk
- www.HidenAnalytical.com

Membrane inlet mass spectrometry method (REOX/MIMS) to measure 15N-nitrate in isotope-enrichment experiments

Authors: Xianbiao Lin^{1,2,3,5‡}, Kaijun Lu^{2‡}, Amber K. Hardison^{2,4}, Zhanfei Liu², Xin Xu², Dengzhou Gao³, Jun Gong^{1,5*}, Wayne S. Gardner²

¹Laboratory of Microbial Ecology and Matter Cycles, School of Marine Sciences, Sun Yat-Sen University, Zhuhai 519082, China

²The University of Texas at Austin Marine Science Institute, 750 Channel View Drive, Port Aransas, Texas 78373, USA

³School of Geographic Sciences, Key Laboratory of Geographic Information Science of the Ministry of Education, East China Normal University, Shanghai, 200241, China

⁴Virginia Institute of Marine Sciences, College of William and Mary, Gloucester Point, Virginia, 23062, USA ⁵Southern Laboratory of Ocean Science and Engineering (Guangdong, Zhuhai), Zhuhai, 519000, China

Using ¹⁵N stable isotope as a tracer to quantify N transformation rates in isotope-enrichment experiments improves understanding of the N cycle in various ecosystems. However, measuring ¹⁵Nnitrate (15NO₃-) in small volumes of water for these experiments is a major challenge due to the inconvenience of preparing samples by traditional techniques. Lin et al. (2021) developed a "REOX/MIMS" method by applying membrane inlet mass spectrometry (MIMS, Hiden HPR-40, Hiden **Analytical Ltd., Warrington, UK)** to determining ¹⁵NO₃⁻ concentrations in a small volume of water from isotope-enrichment experiments after converting the dissolved inorganic N to N₂. The nitrates $(NO_3^- + NO_2^-)$ were reduced to NH_4^+ with zinc powder, and the ammonium (NH_4^+) was then oxidized to N_2 by hypobromite iodine solution. The resulting $^{29}N_2$ and $^{30}N_2$ were measured via MIMS. This optimized protocol provides a sensitive (\sim 0.1 μ M) and precise (relative standard deviation = 0.1– 4.37%) approach to quantify $^{15}NO_3$ concentrations (0.1–500 μ M) in water samples over a wide range of salinities (0–35 ‰) and in 2 M KCl solution with excellent calibration curves ($R^2 \ge 0.9996$, p < 0.99960.0001). The method was combined with ¹⁵NO₃⁻ isotope-enrichment incubation experiments to measure gross nitrification and gross NO₃⁻ immobilization rates in various ecosystems. It was rapid, accurate, and cost-effective. Future applications of this efficient approach will inform scientists, modelers and decision makers about mechanisms, sources, fates, and effects of NO₃⁻ delivered to or produced in numerous aquatic and terrestrial ecosystems. The privous OX/MIMS method has be expanded successfully to determine DNRA, N fixation, mineralization and immobilization with isotope tracer or dilution techniques in sediments of aquatic environments (Yin, et al, 2014; Hou et al., 2018; Lin et al., 2016a; Lin et al., 2016b; Lin et al., 2017). Extending use of the MIMS to measurements of $^{15}\text{NO}_3^-$ using the described REOX/MIMS method provides a convenient and cost-effective approach to determine gross nitrification, gross N mineralization and immobilization, and N fixation rates, and extends the application field of MIMS to terrestrial ecosystems. In addition, REOX/MIMS has been extended to determine DO15N concentration using UV oxidation (Lu et al., 2020). With this method, we can quantify the main N-transformation processes in the soils/sediments from various ecosystems using a MIMS.






Fig. 1. The general procedure of the "REOX/MIMS" method for determination of $^{15}NO_3^-$ in aqueous samples (a) and the schematic diagram of self-assembled membrane injection mass spectrometry system (b); The main components of this system are: sample vial (S), injection peristaltic pump (P), constant temperature water bath (T), stainless steel capillary (SC), membrane injector (M, including a gas-permeable silicone elastomer tube and a thick glass tubing), waste recovery bottle (W), cold trap (C), copper reduction furnace (F, containing a quartz tube with reduced copper wire), vacuum system (V), Ion source (I), quadrupole mass analyzer (Q), secondary electron multiplier (SEM), data processing system (D).




Fig. 2. Relationships of the known $^{15}NO_3^-$ concentrations with measured signal intensities of total ^{15}N ($^{29}N_2+2\times^{30}N_2$) under optimal condition at salinity of 0, 15, and 35 ‰, as well as at solution of 2 M KCl. Vertical bars denote the standard errors (n = 3). S and P represent salinity and 2 M KCl, respectively.






Fig. 3. (a) Relationships of the measured ¹⁵N fraction (0.5, 1, 5, 10, 25, 50, and 75%) with expected ¹⁵N fraction for standards at 500 μ M, (b) ¹⁵N abundances measured by at different NO₃⁻ concentrations (n = 3, mean and standard deviation)

References

- 1) Lin XB, Lu KJ, Hardison AK, Xu X, Gao DZ, Gong J*, Gardner WS. New membrane inlet mass spectrometry method (REOX/MIMS) to measure ¹⁵N-nitrate in isotope-enrichment experiments: Analytical technique and applications. Ecological Indicartors, 2021, 126: 107639.
- 2) Lin XB, Hou LJ, Liu M, Li XF, Yin GY, Zheng YL, et al. Gross nitrogen mineralization in surface sediments of the Yangtze Estuary. PLoS ONE 2016a; 11: e0151930. Doi: 10.1371/journal.pone.0151930.
- 3) Lin XB, Hou LJ, Liu M, Li XF, Zheng YL, Yin GY, et al. Nitrogen mineralization and immobilization in sediments of the East China Sea: Spatiotemporal variations and environmental implications. J. Geophys. Res. Biogeosci. 2016b; 121: 2842-2855. Doi: 10.1002/2016JG003499.
- 4) Lin XB, Li XF, Gao DZ, Liu M, Cheng L. Ammonium production and removal in the sediments of Shanghai River networks: Spatiotemporal variations, controlling factors, andervironmental implications. J. Geophys. Res. Biogeosci. 2017a; 122: 2461-2478. Doi: 10.1002/2017JG003769.
- 5) Yin GY, Hou LJ, Liu M, Liu ZF, Gardner WS. A novel membrane inlet mass spectrometer method to measure 15 NH₄ $^+$ for isotope-enrichment experiments in aquatic ecosystems. Environ. Sci. Technol. 2014; 48: 9555-9562. Doi: 10.1021/es501261s.
- 6) Lu KJ, Lin XB, Gardner WS, Liu ZF. A streamlined method to quantify the fates of ¹⁵N in seawater samples amended with ¹⁵N-labeled organic nitrogen. Limnology and Oceanography: Methods. 2020; 18(2): 52-62. Doi:10.1002/lom3.10345.

Paper Reference: "Membrane inlet mass spectrometry method (REOX/MIMS) to measure ¹⁵N-nitrate in isotope-enrichment experiments", *Ecological Indicators* **2021**, 126, 107639 https://doi.org/10.1016/j.ecolind.2021.107639

Hiden Product: HPR-40 DSA