

New Literature Information



New Hidden Product Catalogue: Catalysis and Thermal Analysis

Hidden Analytical release a new catalogue detailing their latest mass spectrometers and micro-reactors for analysis of gaseous thermal processes. The series includes systems for application to process pressures from 5 bar through to ultrahigh vacuum (UHV).

The CATLAB 1000°C micro-reactors with integrated mass spectrometer enable detailed quantification of catalytic and thermal process activity with full control of temperature and temperature ramp, feed gas species selection and flow rate. Additional process data such as sample mass may also be integrated. The TPD Workstation for UHV thermal desorption measurement features a sample heating platform and fast-response high-sensitivity quadrupole mass spectrometer to maximise data accumulation rate, together with a load lock with multiple-sample transfer mechanism to maintain optimum vacuum integrity. The HPR-60 MBMS system with multiple-stage pressure reduction enables measurement of reactive ions and radicals directly from plasma and flame processes at pressures up to 5 bar.

An extensive range of process interfaces enables direct coupling of the Hidden mass spectrometers with user processes including thermal and thermal-gravimetric systems and vacuum furnaces. The unique SPACI-MS interface is developed specifically for in situ analysis within a catalytic monolith for precise performance evaluation.

New Catalysis & Thermal Analysis Catalogue

The catalogue is available to download at <http://tinyurl.com/cta-catalogue-166-3> or you can request your free print copy via our website: <http://tinyurl.com/request-print-copy>. For further information on all Hidden Analytical products contact Hidden Analytical at info@hidden.co.uk or visit the main website at: www.HiddenAnalytical.com.

--- ends ---